

兔抗 TRPV3 多克隆抗体

中文名称：兔抗 TRPV3 多克隆抗体

英文名称：Anti-TRPV3 rabbit ployclonal antibody

别名：OLMS; VRL3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TRPV3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene product belongs to a family of nonselective cation channels that function in a variety of processes, including temperature sensation and vasoregulation. The thermosensitive members of this family are expressed in subsets of sensory neurons that terminate in the skin, and are activated at distinct physiological temperatures. This channel is activated at temperatures between 22 and 40 degrees C. This gene lies in close proximity to another family member gene on chromosome 17, and the two encoded proteins are thought to associate with each other to form heteromeric

	ric channels. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TRPV3
Immunogen:	Synthetic peptide of human TRPV3
Full name:	transient receptor potential cation channel, subfamily V, member 3
Synonyms :	OLMS; VRL3
SwissProt:	Q8NET8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	25-100

