

## 兔抗 RABL2A 多克隆抗体

中文名称: 兔抗 RABL2A 多克隆抗体

英文名称: Anti-RABL2A rabbit polyclonal antibody

储存: 冷冻 (-20℃) 避光

抗原: RABL2A

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	This gene is a member of the RAB gene family which belongs to the RAS GTPase superfamily. The proteins in the family of RAS-related signaling molecules are small GTP-binding proteins that play important roles in the regulation of exocytotic and endocytotic pathways. This gene maps to the site of an ancestral telomere fusion event and may be a subtelomeric gene. Two alternatively spliced transcript variants have been identified, both encoding the same protein.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	RABL2A
<b>Immunogen:</b>	Synthesized peptide derived from internal of hum

	an RABL2A.
<b>Full name:</b>	RAB, member of RAS oncogene family-like 2A
<b>Synonyms :</b>	0
<b>SwissProt:</b>	Q9UBK7
<b>WB Predicted band size:</b>	26 kDa
<b>WB Positive control:</b>	Huvec, 293 and COLO205 cells lysates
<b>WB Recommended dilution:</b>	500-3000

