

兔抗 HTR1B 多克隆抗体

中文名称： 兔抗 HTR1B 多克隆抗体

英文名称： Anti-HTR1B rabbit polyclonal antibody

别名： S12; 5-HT1B; HTR1D2; HTR1DB; 5-HT1DB

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： HTR1B

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The neurotransmitter serotonin (5-hydroxytryptamine; 5-HT) exerts a wide variety of physiologic functions through a multiplicity of receptors and may be involved in human neuropsychiatric disorders such as anxiety, depression, or migraine. These receptors consist of several main groups subdivided into several distinct subtypes on the basis of their pharmacologic characteristics, coupling to intracellular second messengers, and distribution within the nervous system (Zifa and Fillion, 1992 [PubMed 1359584]). The serotonergic receptors belong to the multigene family of receptors coupled to guanine nucleotide-binding proteins.

Applications:	ELISA, IHC
Name of antibody:	HTR1B
Immunogen:	Synthetic peptide of human HTR1B
Full name:	5-hydroxytryptamine (serotonin) receptor 1B, G protein-coupled
Synonyms :	S12; 5-HT1B; HTR1D2; HTR1DB; 5-HT1DB
SwissProt:	P28222
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human tonsil
IHC Recommend dilution:	25-100

