

兔抗 RABL2B 多克隆抗体

中文名称：兔抗 RABL2B 多克隆抗体

英文名称： Anti-RABL2B rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： RABL2B

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The RABL2B protein is a member of the RAB gene family which belongs to the RAS GTPase superfamily. RABL2B is located within a subtelomeric region of 22q13.3. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	RABL2B
Immunogen:	Fusion protein of human RABL2B
Full name:	RAB, member of RAS oncogene family-like 2B
SwissProt:	Q9UNT1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human breast cancer

IHC Recommend dilution:

25-100

