

兔抗 NCF1(Ab-304) 多克隆抗体

中文名称: 兔抗 NCF1(Ab-304) 多克隆抗体

英文名称: Anti-NCF1(Ab-304) rabbit polyclonal antibody

别名: NCF1A; NOXO2; p47phox; SH3PXD1A

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

抗原: NCF1(Ab-304)

反应种属: Human


标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	The protein encoded by this gene is a 47 kDa cytosolic subunit of neutrophil NADPH oxidase. This oxidase is a multicomponent enzyme that is activated to produce superoxide anion. Mutations in this gene have been associated with chronic granulomatous disease.
Applications:	WB, IF
Name of antibody:	NCF1(Ab-304)
Immunogen:	Synthesized non-phosphopeptide derived from human Neutrophil Cytosol Factor 1 around the phosphorylation site of serine 304 (R-S-S(p)-I-R).
Full name:	neutrophil cytosolic factor 1

Synonyms :	NCF1A; NOXO2; p47phox; SH3PXD1A
SwissProt:	P14598
WB Predicted band size:	45 kDa
WB Positive control:	COLO cell lysate
WB Recommended dilution:	500-3000
IF positive control:	Hela cells
IF Recommend dilution:	100-500



mlbio 梅联生物
Good elisakit producers

