

兔抗 TSEN2 多克隆抗体

中文名称:	兔抗 TSEN2 多克隆抗体
英文名称:	Anti-TSEN2 rabbit polyclonal antibody
别名:	tRNA splicing endonuclease subunit 2; SEN2; PCH2B; SEN2L
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	TSEN2
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	This gene encodes one of the subunits of the tRNA splicing endonuclease. This endonuclease catalyzes the first step in RNA splicing which is the removal of introns. Mutations in this gene have been associated with pontocerebellar hypoplasia type 2. A pseudogene has been identified on chromosome 4. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TSEN2
Immunogen:	Fusion protein of human TSEN2
Full name:	tRNA splicing endonuclease subunit 2

Synonyms:	SEN2; PCH2B; SEN2L
SwissProt:	Q8NCEO
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	50-300

