

兔抗 NCF1(Phospho-Ser359) 多克隆抗体

中文名称：兔抗 NCF1(Phospho-Ser359) 多克隆抗体

英文名称：Anti-NCF1(Phospho-Ser359) rabbit polyclonal antibody

别名：NCF1A; NOXO2; p47phox; SH3PXD1A

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：NCF1(Phospho-Ser359)

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The protein encoded by this gene is a 47 kDa cytosolic sub unit of neutrophil NADPH oxidase. This oxidase is a multicomponent enzyme that is activated to produce superoxide anion. Mutations in this gene have been associated with chronic granulomatous disease.
Applications:	WB
Name of antibody:	NCF1(Phospho-Ser359)
Immunogen:	Peptide sequence around phosphorylation site of Serine 359 (Q-R-S(p)-K-P) derived from Human p47 phox.
Full name:	neutrophil cytosolic factor 1
Synonyms :	NCF1A; NOXO2; p47phox; SH3PXD1A

SwissProt:	P14598
WB Predicted band size:	45 kDa
WB Positive control:	Hela cells lysates
WB Recommended dilution:	500-1000

