

## 兔抗 GST Tag 多克隆抗体

- 中文名称： 兔抗 GST Tag 多克隆抗体
- 英文名称： Anti-GST Tag rabbit polyclonal antibody
- 相关类别： 标签抗体
- 储 存： 冷冻（-20℃）
- 宿 主： Rabbit
- 抗 原： GST Tag
- 反应种属： Fusion proteins containing the GST Tag
- 标 记 物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	GST Tag
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	Glutathione S-Transferase
<b>ELISA Recommended dilution:</b>	1000-10000
<b>WB Predicted band size:</b>	63 kDa

<b>WB Positive control:</b>	Recombinant GST tag protein
<b>WB Recommended dilution:</b>	1000-8000

