

兔抗 TSEN15 多克隆抗体

中文名称: 兔抗 TSEN15 多克隆抗体

英文名称: Anti-TSEN15 rabbit polyclonal antibody

别名: sen15; C1orf19

抗原: TSEN15

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a subunit of the tRNA splicing endonuclease, which catalyzes the removal of introns from tRNA precursors. Alternative splicing results in multiple transcript variants. There is a pseudogene of this gene on chromosome 17.
Applications:	ELISA, IHC
Name of antibody:	TSEN15
Immunogen:	Full length fusion protein
Full name:	tRNA splicing endonuclease subunit 15
Synonyms:	sen15; C1orf19
SwissProt:	Q8WW01

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and human breast cancer
IHC Recommend dilution:	25-100



