

## 兔抗 NCK1 多克隆抗体

中文名称：兔抗 NCK1 多克隆抗体

英文名称： Anti-NCK1 rabbit polyclonal antibody

别名： NCK, nck-1, NCKalpha

抗原： NCK1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is one of the signaling and transforming proteins containing Src homology 2 and 3 (SH2 and SH3) domains. It is located in the cytoplasm and is an adaptor protein involved in transducing signals from receptor tyrosine kinases to downstream signal recipients such as RAS. Alternatively spliced transcript variants encoding different isoforms have been found.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NCK1
<b>Immunogen:</b>	Fusion protein of human NCK1
<b>Full name:</b>	NCK adaptor protein 1

<b>Synonyms :</b>	NCK, nck-1, NCKalpha
<b>SwissProt:</b>	P16333
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human lung cancer and human colon cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	43 kDa
<b>WB Positive control:</b>	Jurkat and K562 cells
<b>WB Recommended dilution:</b>	200-1000



