

兔抗 HTR4 多克隆抗体

中文名称：兔抗 HTR4 多克隆抗体

英文名称：Anti-HTR4 rabbit polyclonal antibody

别名：5-HT4, 5-HT4R

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HTR4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the family of serotonin receptors, which are G protein coupled receptors that stimulate cAMP production in response to serotonin (5-hydroxytryptamine). The gene product is a glycosylated transmembrane protein that functions in both the peripheral and central nervous system to modulate the release of various neurotransmitters. Multiple transcript variants encoding proteins with distinct C-terminal sequences have been described.
Applications:	ELISA, WB
Name of antibody:	HTR4

Immunogen:	Synthetic peptide of humanHTR4
Full name:	5-hydroxytryptamine (serotonin) receptor 4, G protein-coupled
Synonyms :	5-HT4; 5-HT4R
SwissProt:	Q13639
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	44 kDa
WB Positive control:	Hela, 293T, A549 and Jurkat cells
WB Recommended dilution:	200-1000

