

兔抗 TSHR 多克隆抗体

中文名称：兔抗 TSHR 多克隆抗体

英文名称： Anti-TSHR rabbit polyclonal antibody

别名： LGR3; CHNG1; hTSHR-I

储存： 冷冻（-20℃）

抗原： TSHR

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a membrane protein and a major controller of thyroid cell metabolism. The encoded protein is a receptor for thyrothropin and thyrostimulin, and its activity is mediated by adenylate cyclase. Defects in this gene are a cause of several types of hyperthyroidism. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TSHR
Immunogen:	Fusion protein of human TSHR
Full name:	thyroid stimulating hormone receptor

Synonyms :	LGR3; CHNG1; hTSHR-I
SwissProt:	P16473
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

