

兔抗 RADIL 多克隆抗体

中文名称：兔抗 RADIL 多克隆抗体

英文名称： Anti-RADIL rabbit polyclonal antibody

别名： Rap associating with DIL domain; RASIP2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RADIL

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Downstream effector of Rap required for cell adhesion and migration of neural crest precursors during development.
Applications:	ELISA, IHC
Name of antibody:	RADIL
Immunogen:	Fusion protein of human RADIL
Full name:	Rap associating with DIL domain
Synonyms:	RASIP2
SwissProt:	Q96JH8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer and Human tonsil

IHC Recommend dilution:

40-200



