

## 小鼠抗 ECFP Tag 单克隆抗体

中文名称: 小鼠抗 ECFP Tag 单克隆抗体

英文名称: Anti-ECFP Tag mouse monoclonal antibody

别名: Enhanced Cyan Fluorescent Protein

相关类别: 标签抗体

储存: 冷冻 (-20℃)

宿主: Mouse

抗原: ECFP Tag

反应种属: Fusion proteins containing ECFP

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jelly fish <i>Aequorea victoria</i> .
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ECFP Tag
<b>Immunogen:</b>	Fusion protein of ECFP
<b>Full name:</b>	Enhanced Cyan Fluorescent Protein
<b>WB Positive control:</b>	Recombinant ECFP tag protein
<b>WB Recommended dilution:</b>	5000-20000

