

## 兔抗 NCOA6 多克隆抗体

中文名称: 兔抗 NCOA6 多克隆抗体

英文名称: Anti-NCOA6 rabbit polyclonal antibody

别名: NRC; AIB3; ASC2; PRIP; TRBP; RAP250

储存: 冷冻 (-20℃)

抗原: NCOA6

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a transcriptional coactivator that can interact with nuclear hormone receptors to enhance their transcriptional activator functions. This protein has been shown to be involved in the hormone-dependent coactivation of several receptors, including prostanoic acid, retinoid, vitamin D3, thyroid hormone, and steroid receptors. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NCOA6
<b>Immunogen:</b>	Synthetic peptide of human NCOA6

<b>Full name:</b>	nuclear receptor coactivator 6
<b>Synonyms:</b>	NRC; AIB3; ASC2; PRIP; TRBP; RAP250
<b>SwissProt:</b>	Q14686
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

