

兔抗 NCOA7 多克隆抗体

中文名称: 兔抗 NCOA7 多克隆抗体

英文名称: Anti-NCOA7 rabbit polyclonal antibody

别名: ESNA1; TLDC4; ERAP140; NCOA7-AS; Nbla00052; Nbla10993; dJ187J11.3

储存: 冷冻 (-20℃) 避光

抗原: NCOA7

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Enhances the transcriptional activities of several nuclear receptors. Involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA.
Applications:	WB
Name of antibody:	NCOA7
Immunogen:	Synthesized peptide derived from internal of human NCOA7.
Full name:	nuclear receptor coactivator 7
Synonyms :	ESNA1; TLDC4; ERAP140; NCOA7-AS; Nbla00052; Nbla10993; dJ187J11.3
SwissProt:	Q8NI08

WB Predicted band size:	106 kDa
WB Positive control:	LOVO cells, NIH/3T3 cells and RAW264.7 cells lysates
WB Recommended dilution:	500-3000

