

兔抗 NCOR1 多克隆抗体

中文名称：兔抗 NCOR1 多克隆抗体

英文名称： Anti-NCOR1 rabbit polyclonal antibody

别名： N-CoR, TRAC1, N-CoR1, hN-CoR

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NCOR1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a protein that mediates ligand-independent transcription repression of thyroid-hormone and retinoic-acid receptors by promoting chromatin condensation and preventing access of the transcription machinery. It is part of a complex which also includes histone deacetylases and transcriptional regulators similar to the yeast protein Sin3p. This gene is located between the Charcot-Marie-Tooth and Smith-Magenis syndrome critical regions on chromosome 17. Alternate splicing results in multiple transcript variants. Pseudogenes of this gene are found on chromosomes 17 and 20.

Applications:	ELISA, IHC
Name of antibody:	NCOR1
Immunogen:	Synthetic peptide of human NCOR1
Full name:	nuclear receptor corepressor 1
Synonyms :	N-CoR, TRAC1, N-CoR1, hN-CoR
SwissProt:	O75376
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and Human tonsil
IHC Recommend dilution:	100-300

