

兔抗 TSPY1/TSPY3 多克隆抗体

中文名称：兔抗 TSPY1/TSPY3 多克隆抗体

英文名称： Anti-TSPY1/TSPY3 rabbit polyclonal antibody

别名： CT78; TSPY; DYS14; pJA923/CT78

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TSPY1/TSPY3

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is found only in testicular tissue and may be involved in spermatogenesis. Approximately 35 copies of this gene are present in humans, but only a single, nonfunctional orthologous gene is found in mouse. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	TSPY1/TSPY3
Immunogen:	Synthetic peptide of human TSPY1/TSPY3
Full name:	testis specific protein, Y-linked 1/testis specific protein, Y-linked 3

Synonyms :	CT78; TSPY; DYS14; pJA923/CT78
SwissProt:	Q01534/P0CV98
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	30-150

