

兔抗 TSPYL1 多克隆抗体

中文名称：兔抗 TSPYL1 多克隆抗体

英文名称：Anti-TSPYL1 rabbit polyclonal antibody

别名：TSPYL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TSPYL1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is found in the nucleolus and is similar to that of a family of genes on the Y-chromosome. This gene is intronless. Defects in this gene are a cause of sudden infant death with dysgenesis of the testes syndrome (SIDDT).
Applications:	ELISA, WB
Name of antibody:	TSPYL1
Immunogen:	Synthetic peptide of humanTSPYL1
Full name:	TSPY-like 1
Synonyms :	TSPYL
SwissProt:	Q9H0U9

ELISA Recommended dilution:	2000-5000
WB Predicted band size:	49 kDa
WB Positive control:	Mouse brain and heart tissue
WB Recommended dilution:	500-2000

