

## 兔抗 TSPYL2 多克隆抗体

中文名称：兔抗 TSPYL2 多克隆抗体

英文名称：Anti-TSPYL2 rabbit polyclonal antibody

别名：CDA1; CTCL; NP79; TSPX; CINAP; DENTT; SE204; HRIHFB2216

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TSPYL2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the testis-specific protein Y-encoded, TSPY-like/SET/nucleosome assembly protein-1 superfamily. The encoded protein is localized to the nucleolus where it functions in chromatin remodeling and as an inhibitor of cell-cycle progression. This protein may play a role in the suppression of tumor growth.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	TSPYL2
<b>Immunogen:</b>	Synthetic peptide of human TSPYL2
<b>Full name:</b>	TSPY-like 2

<b>Synonyms :</b>	CDA1; CTCL; NP79; TSPX; CINAP; DENTT; SE204; HRIHFB 2216
<b>SwissProt:</b>	Q9H2G4
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	231 cells
<b>WB Recommended dilution:</b>	500-2000

