

## 兔抗 HTR6 多克隆抗体

中文名称: 兔抗 HTR6 多克隆抗体

英文名称: Anti-HTR6 rabbit polyclonal antibody

别名: 5-HT6; 5-HT6R

抗原: HTR6

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a protein that belongs to the seven-transmembrane G protein-coupled receptor family of proteins. The encoded protein couples with the Gs alpha subunit and stimulates a denylate cyclase to activate the cyclic AMP-dependent signaling pathway. This receptor is thought to regulate cholinergic neuronal transmission in the brain. Several antidepressants and antipsychotic drugs have a high affinity for this receptor.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	HTR6
<b>Immunogen:</b>	Synthesized peptide derived from internal of human 5-HT-6.

<b>Full name:</b>	5-hydroxytryptamine (serotonin) receptor 6, G protein-coupled
<b>Synonyms :</b>	5-HT6; 5-HT6R
<b>SwissProt:</b>	P50406
<b>WB Predicted band size :</b>	47 kDa
<b>WB Positive control:</b>	A549 cells lysate
<b>WB Recommended dilution:</b>	500-3000

