

兔抗 NCR3 多克隆抗体

中文名称：兔抗 NCR3 多克隆抗体

英文名称：Anti-NCR3 rabbit polyclonal antibody

别名：1C7; MALS; CD337; LY117; NKp30

相关类别：一抗

抗原：NCR3

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a natural cytotoxicity receptor (NCR) that may aid NK cells in the lysis of tumor cells. The encoded protein interacts with CD3-zeta (CD247), a T-cell receptor. A single nucleotide polymorphism in the 5' untranslated region of this gene has been associated with mild malaria susceptibility. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	NCR3
Immunogen:	Fusion protein of human NCR3

Full name:	natural cytotoxicity triggering receptor 3
Synonyms :	1C7; MALS; CD337; LY117; NKp30
SwissProt:	O14931
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

