

兔抗 NDEL1 多克隆抗体

中文名称：兔抗 NDEL1 多克隆抗体

英文名称：Anti-NDEL1 rabbit polyclonal antibody

别名：EOPA; NDE2; NUDEL; MITAP1; NDE1L1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NDEL1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a coiled-coil protein that plays a role in multiple processes including cytoskeletal organization, cell signaling and neuron migration, outgrowth and maintenance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome X.
Applications:	ELISA, IHC
Name of antibody:	NDEL1
Immunogen:	Fusion protein of human NDEL1
Full name:	nudE neurodevelopment protein 1-like 1

Synonyms :	EOPA; NDE2; NUDEL; MITAP1; NDE1L1
SwissProt:	Q9GZM8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	50-200

