

兔抗 TSTA3 多克隆抗体

中文名称：兔抗 TSTA3 多克隆抗体

英文名称：Anti-TSTA3 rabbit polyclonal antibody

别名：tissue specific transplantation antigen P35B; FX; P35B; SDR4E1

储存：冷冻（-20℃）

抗原：TSTA3

宿主：Rabbit

反应种属：Human, Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Tissue specific transplantation antigen P35B is a NADP(H)-binding protein. It catalyze the two-step epimerase and the reductase reactions in GDP-D-mannose metabolism, converting GDP-4-keto-6-D-deoxymannose to GDP-L-fucose. GDP-L-fucose is the substrate of several fucosyltransferases involved in the expression of many glycoconjugates, including blood group ABH antigens and developmental adhesion antigens. Mutations in this gene may cause leukocyte adhesion deficiency, type II.
Applications:	ELISA, IHC
Name of antibody:	TSTA3
Immunogen:	Fusion protein of human TSTA3
Full name:	tissue specific transplantation antigen P35B

Synonyms:	FX; P35B; SDR4E1
SwissProt:	Q13630
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human brain
IHC Recommend dilution:	50-100



