

## 兔抗 NDFIP1 多克隆抗体

中文名称：兔抗 NDFIP1 多克隆抗体

英文名称：Anti-NDFIP1 rabbit polyclonal antibody

别名：N4WBP5

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：NDFIP1

反应种属：Human Mouse Rat

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	N4WBP5
<b>Applications:</b>	WB
<b>Name of antibody:</b>	NDFIP1
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human NDFIP1.
<b>Full name:</b>	Nedd4 family interacting protein 1
<b>Synonyms :</b>	N4WBP5
<b>SwissProt:</b>	Q9BT67
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	293 cells lysate
<b>WB Recommended dilution:</b>	500-3000

