

兔抗 RAP1A 多克隆抗体

中文名称：兔抗 RAP1A 多克隆抗体

英文名称： Anti-RAP1A rabbit polyclonal antibody

别名： RAP1; KREV1; KREV-1; SMGP21

抗原： RAP1A

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the Ras family of small GTPases. The encoded protein undergoes a change in conformational state and activity, depending on whether it is bound to GTP or GDP. This protein is activated by several types of guanine nucleotide exchange factors (GEFs), and inactivated by two groups of GTPase-activating proteins (GAPs). The activation status of the encoded protein is therefore affected by the balance of intracellular levels of GEFs and GAPs.
Applications:	ELISA, WB, IHC

Name of antibody:	RAP1A
Immunogen:	Full length fusion protein
Full name:	RAP1A, member of RAS oncogene family
Synonyms :	RAP1; C21KG; G-22K; KREV1; KREV-1; SMGP21
SwissProt:	P62834
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	21 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	200-1000



