

兔抗 NDRG1 多克隆抗体

中文名称：兔抗 NDRG1 多克隆抗体

英文名称：Anti-NDRG1 rabbit polyclonal antibody

别名：GC4, RTP, DRG1, NDR1, NMSL, TDD5, CAP43, CMT4D, DRG-1, HMSNL, RIT42, TARG1, PROXY1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NDRG1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene is a member of the N-myc downregulated gene family which belongs to the alpha/beta hydrolase superfamily. The protein encoded by this gene is a cytoplasmic protein involved in stress responses, hormone responses, cell growth, and differentiation. The encoded protein is necessary for p53-mediated caspase activation and apoptosis. Mutations in this gene are a cause of Charcot-Marie-Tooth disease type 4D, and expression of this gene may be a prognostic indicator for several types of cancer. Alternatively spliced transcr

	ipt variants encoding multiple isoforms have been observed for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	NDRG1
Immunogen:	Synthetic peptide of human NDRG1
Full name:	N-myc downstream regulated 1
Synonyms :	GC4, RTP, DRG1, NDR1, NMSL, TDD5, CAP43, CMT4D, DRG-1, HMSNL, RIT42, TARG1, PROXY1
SwissProt:	Q92597
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain and Human esophagus cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	43 kDa
WB Positive control:	Human lymphoma tissue
WB Recommended dilution:	500-2000



