

## 小鼠抗 GST Tag 单克隆抗体

- 中文名称： 小鼠抗 GST Tag 单克隆抗体
- 英文名称： Anti-GST Tag mouse monoclonal antibody
- 相关类别： 标签抗体
- 储 存： 冷冻（-20℃）
- 宿 主： Mouse
- 抗 原： GST Tag
- 反应种属： Fusion proteins containing the GST Tag
- 标 记 物： Unconjugate
- 克隆类型： mouse monoclonal

### 技术规格

|                                |  |
|--------------------------------|--|
| <b>Background:</b>             | GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut. |
| <b>Applications:</b>           | WB   |
| <b>Name of antibody:</b>       | GST Tag  |
| <b>Immunogen:</b>              | Full length fusion protein   |
| <b>Full name:</b>              | Glutathione S-Transferase  |
| <b>WB Predicted band size:</b> | 26 kDa   |
| <b>WB Positive control:</b>    | Recombinant GST tag protein  |

WB Recommended dilution: 5000-10000

