

小鼠抗 GST 单克隆抗体

中文名称： 小鼠抗 GST 单克隆抗体

英文名称： Anti-GST mouse monoclonal antibody

相关类别： 标签抗体

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： GST

反应种属： GST tag

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.
Applications:	WB
Name of antibody:	GST
Immunogen:	Recombinant GST protein
Full name:	Glutathione S-Transferase
WB Positive control:	six different HEK293T over-expression lysates containing different GST tagged recombinant proteins, and one



	GST-DDK bacterial lysate produced in E.coli.
WB Recommended dilution:	2000-10000