

兔抗 TTC4 多克隆抗体

中文名称: 兔抗 TTC4 多克隆抗体

英文名称: Anti-TTC4 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: TTC4

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains tetratricopeptide (TPR) repeats, which often mediate protein-protein interactions and chaperone activity. The encoded protein interacts with heat shock proteins 70 and 90. Alternative splicing results in multiple transcript variants. Naturally-occurring readthrough transcription occurs from upstream gene MROH (maestro heat-like repeat family member 7) to this gene.
Applications:	ELISA, IHC
Name of antibody:	TTC4
Immunogen:	Fusion protein of human TTC4
Full name:	tetratricopeptide repeat domain 4
SwissProt:	O95801

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer;
IHC Recommend dilution:	40-200

