

小鼠抗 HA tag 单克隆抗体

中文名称： 小鼠抗 HA tag 单克隆抗体

英文名称： Anti-HA tag mouse monoclonal antibody

别名： Human influenza hemagglutinin Tag; HA epitope tag; Hemagglutinin Tag

相关类别： 标签抗体

储存： 冷冻（-20℃）

宿主： Mouse

抗原： HA tag

反应种属： Fusion proteins containing the HA Tag

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Human influenza hemagglutinin (HA) is a surface glyco protein required for the infectivity of the human virus. The HA tag is derived from the HA-molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein. This tag facilitates the detection, isolation, and purification of the proteins.
Applications:	WB

Name of antibody:	HA tag
Immunogen:	Synthetic peptide conjugated with KLH
Full name:	Human influenza hemagglutinin Tag
Synonyms:	HA epitope tag, Hemagglutinin Tag
WB Predicted band size:	55 kDa
WB Positive control:	Fusion proteins containing the HA Tag
WB Recommended dilution:	10000-50000

