

## 兔抗 TTC12 多克隆抗体

中文名称: 兔抗 TTC12 多克隆抗体

英文名称: Anti-TTC12 rabbit polyclonal antibody

别名: tetratricopeptide repeat domain 12; TPARM

抗原: TTC12

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate


克隆类型: rabbit polyclonal

### 技术规格

**Background:**

TTC12 (tetratricopeptide repeat domain 12), also known as TPARM, is a 705 amino acid cytoplasmic protein containing a TPR repeat domain and three armadillo repeat motifs. TTC12 is expressed in testis, prostate, lung, germinal center B cells, neuroblastoma, teratocarcinoma, colon cancer and gastric cancer. The gene encoding TTC12 is located in a region of human chromosome 11 that is commonly deleted in a variety of cancers, therefore, it is considered a candidate tumor suppressor gene. Human chromosome 11 houses over 1,400 genes and comprises nearly 4% of the human genome.

<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	TTC12
<b>Immunogen:</b>	Fusion protein of human TTC12
<b>Full name:</b>	tetratricopeptide repeat domain 12
<b>Synonyms:</b>	TPARM
<b>SwissProt:</b>	Q9H892
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	Human bladder transitional cell carcinoma grade 2-3 tissue lysate
<b>WB Recommended dilution:</b>	500-2000



**mlbio** 梅联生物  
Good elisakit producers

