

兔抗 TTC16 多克隆抗体

中文名称：兔抗 TTC16 多克隆抗体

英文名称：Anti-TTC16 rabbit polyclonal antibody

别名：tetratricopeptide repeat domain 16

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TTC16

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, WB
Name of antibody:	TTC16
Immunogen:	Fusion protein of human TTC16
Full name:	tetratricopeptide repeat domain 16
SwissProt:	Q8NEE8
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	98 kDa
WB Positive control:	231 cell lysate
WB Recommended dilution:	1000-5000

