

## 兔抗 TTF2 多克隆抗体

中文名称：兔抗 TTF2 多克隆抗体

英文名称： Anti-TTF2 rabbit polyclonal antibody

别名： HuF2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TTF2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the SWI2/SNF2 family of proteins, which play a critical role in altering protein-DNA interactions. The encoded protein has been shown to have dsDNA-dependent ATPase activity and RNA polymerase II termination activity. This protein interacts with cell division cycle 5-like, associates with human splicing complexes, and plays a role in pre-mRNA splicing.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TTF2
<b>Immunogen:</b>	Synthetic peptide of human TTF2
<b>Full name:</b>	transcription termination factor, RNA polymerase II

<b>Synonyms :</b>	HuF2
<b>SwissProt:</b>	Q9UNY4
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	5-20

