



兔抗 HUWE1 多克隆抗体

中文名称： 兔抗 HUWE1 多克隆抗体

英文名称： Anti-HUWE1 rabbit polyclonal antibody

级 别： MULE; Ib772; LASU1; UREB1; HECTH9; URE-B1; ARF-BP1; HSPC272

抗 原： HUWE1

储 存： 冷冻 (-20°C)

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat


标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein containing a C-terminal HECT (E6AP type E3 ubiquitin protein ligase) domain that functions as an E3 ubiquitin ligase. The encoded protein is required for the ubiquitination and subsequent degradation of the anti-apoptotic protein Mcl1 (myeloid cell leukemia sequence 1 (BCL2-related)). This protein also ubiquitinates the p53 tumor suppressor, core histones, and DNA polymerase beta. Mutations in this gene are associated with Turner type X-linked syndromic mental retardation.
Applications:	ELISA, IHC
Name of antibody:	HUWE1

Immunogen:	Synthetic peptide of human HUWE1
Full name:	HECT, UBA and WWE domain containing 1, E3 ubiquitin protein ligase
Synonyms:	MULE; Ib772; LASU1; UREB1; HECTH9; URE-B1; ARF-BP1; HSPC272
SwissProt:	Q7Z6Z7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and human lung cancer
IHC Recommend dilution:	10-50



mlbio 瑪聯生物
Good elisakit producers

