

兔抗 TTI1 多克隆抗体

中文名称：兔抗 TTI1 多克隆抗体

英文名称： Anti-TTI1 rabbit polyclonal antibody

别名： smg-10; KIAA0406

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TTI1

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

Regulator of the DNA damage response (DDR). Part of the TTT complex that is required to stabilize protein levels of the phosphatidylinositol 3-kinase-related protein kinase (PIKK) family proteins. The TTT complex is involved in the cellular resistance to DNA damage stresses, like ionizing radiation (IR), ultraviolet (UV) and mitomycin C (MMC). Together with the TTT complex and HSP90 may participate in the proper folding of newly synthesized PIKKs. Promotes assembly, stabilizes and maintains the activity of mTORC1 and mTORC2 complexes, which regulate cell growth and survival in re

	response to nutrient and hormonal signals.
Applications:	ELISA, IHC
Name of antibody:	TTI1
Immunogen:	Synthetic peptide of humanTTI1
Full name:	TELO2 interacting protein 1
Synonyms :	smg-10; KIAA0406
SwissProt:	O43156
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human prostate cancer and Human esophagus cancer
IHC Recommend dilution:	30-150

