

小鼠抗 mStrawberry 单克隆抗体

中文名称： 小鼠抗 mStrawberry 单克隆抗体

英文名称： Anti-mStrawberry mouse monoclonal antibody

相关类别： 标签抗体

储 存： 冷冻（-20℃） 避光

宿 主： Mouse

抗 原： mStrawberry

反应种属： mStrawberry

标 记 物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of Mrfp. The mStrawberry fluorescent protein photobleaches rapidly.
Applications:	WB
Name of antibody:	mStrawberry
Immunogen:	Purified mStrawberry protein expressed in HEK293 cell
WB Predicted band size:	27 kDa
WB Positive control:	HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY mStrawberry ne) cDNA for 48 hrs and lysed
WB Recommended dilution:	1000-5000