

兔抗 TTLL1 多克隆抗体

中文名称：兔抗 TTLL1 多克隆抗体

英文名称：Anti-TTLL1 rabbit polyclonal antibody

别名：tubulin tyrosine ligase like 1; C22orf7; HS323M22B

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TTLL1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Catalytic subunit of the neuronal tubulin polyglutamylase complex. Modifies alpha- and beta-tubulin, generating side chains of glutamate on the gamma-carboxyl groups of specific glutamate residues within the C-terminal tail of alpha- and beta-tubulin (By similarity).
Applications:	ELISA, IHC
Name of antibody:	TTLL1
Immunogen:	Fusion protein of human TTLL1
Full name:	tubulin tyrosine ligase like 1
Synonyms:	C22orf7; HS323M22B
SwissProt:	O95922

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

