

兔抗 NDUFA4L2 多克隆抗体

中文名称: 兔抗 NDUFA4L2 多克隆抗体

英文名称: Anti-NDUFA4L2 rabbit polyclonal antibody

别名: NUOMS

储存: 冷冻 (-20℃) 避光

抗原: NDUFA4L2

宿主: Rabbit

相关类别: 一抗

反应种属: Human Mouse

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

| | |
|--------------------------------|---|
| Background: | Belongs to the complex I NDUFA4 subunit family. |
| Applications: | WB |
| Name of antibody: | NDUFA4L2 |
| Immunogen: | Synthesized peptide derived from internal of human NDUFA4L2. |
| Full name: | NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 4-like 2 |
| Synonyms : | NUOMS |
| SwissProt: | Q9NRX3 |
| WB Predicted band size: | 10 kDa |
| WB Positive control: | Jurkat cells, COLO cells, Huvec cells and HepG2 cells lysates |

WB Recommended dilution: 500-3000

