

兔抗 RASGEF1A 多克隆抗体

中文名称：兔抗 RASGEF1A 多克隆抗体

英文名称： Anti-RASGEF1A rabbit polyclonal antibody

别名： CG4853

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RASGEF1A

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

RASGEF1A (RasGEF domain family, member 1A) is a 481 amino acid guanine-nucleotide exchange factor (GEF) for H-Ras, K-Ras, and N-Ras. RASGEF1A may also play a role in cell migration. RASFE1A is highly expressed in spinal chord and brain and contains one N-terminal Ras-GEF domain and one Ras-GEF domain. The gene that encodes RASGEF1A maps to human chromosome 10, which contains over 800 genes and 135 million nucleotides, making up nearly 4.5% of the human genome. PTEN is an important tumor suppressor gene located on chromosome 10 and, w

	hen defective, causes a genetic predisposition to cancer development known as Cowden syndrome.
Applications:	ELISA, IHC
Name of antibody:	RASGEF1A
Immunogen:	Fusion protein of human RASGEF1A
Full name:	RasGEF domain family, member 1A
Synonyms :	CG4853
SwissProt:	Q8N9B8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human thyroid cancer
IHC Recommend dilution:	20-100

