

兔抗 IBA57 多克隆抗体

中文名称：兔抗 IBA57 多克隆抗体

英文名称：Anti-IBA57 rabbit polyclonal antibody

别名：MMDS3; SPG74; C1orf69

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IBA57

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene localizes to the mitochondrion and is part of the iron-sulfur cluster assembly pathway. The encoded protein functions late in the biosynthesis of mitochondrial 4Fe-4S proteins. Defects in this gene have been associated with autosomal recessive spastic paraplegia-74 and with multiple mitochondrial dysfunctions syndrome-3. Two transcript variants encoding different isoforms have been found for this gene. The smaller isoform is not likely to be localized to the mitochondrion since it lacks the amino-terminal transit peptide.
Applications:	ELISA, IHC
Name of antibody:	IBA57

Immunogen:	Synthetic peptide of human IBA57
Full name:	IBA57 homolog, iron-sulfur cluster assembly
Synonyms:	MMDS3; SPG74; C1orf69
SwissProt:	Q5T440
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and human brain
IHC Recommend dilution:	30-150



