

兔抗 TUBA1B 多克隆抗体

中文名称：兔抗 TUBA1B 多克隆抗体

英文名称：Anti-TUBA1B rabbit polyclonal antibody

别名：tubulin alpha 1b; K-ALPHA-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TUBA1B

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.
Applications:	ELISA, WB, IHC
Name of antibody:	TUBA1B
Immunogen:	Fusion protein of human TUBA1B
Full name:	tubulin alpha 1b
Synonyms:	K-ALPHA-1
SwissProt:	P68363
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-300
WB Predicted band size:	50 kDa
WB Positive control:	NIH/3T3 cell lysate
WB Recommended dilution:	200-1000



