

兔抗 RASSF6 多克隆抗体

中文名称：兔抗 RASSF6 多克隆抗体

英文名称： Anti-RASSF6 rabbit polyclonal antibody

储 存： 冷冻（-20℃） 避光

抗 原： RASSF6

宿 主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the Ras-association domain family (RASSF). Members of this family form the core of a highly conserved tumor suppressor network, the Salvador-Warts-Hippo (SWH) pathway. The protein encoded by this gene is a Ras effector protein that induces apoptosis. A genomic region containing this gene has been linked to susceptibility to viral bronchiolitis. Alternative splicing results in multiple transcript variants and protein isoforms.
Applications:	WB
Name of antibody:	RASSF6
Immunogen:	Synthesized peptide derived from internal of human RASSF6.

Full name:	Ras association (RalGDS/AF-6) domain family member 6
Synonyms :	0
SwissProt:	Q6ZTQ3
WB Predicted band size:	43 kDa
WB Positive control:	Hela cells lysate
WB Recommended dilution:	500-3000

