

## 兔抗 ICAM2 多克隆抗体

中文名称：兔抗 ICAM2 多克隆抗体

英文名称：Anti-ICAM2 rabbit polyclonal antibody

抗 原：ICAM2

别 名：CD102

储 存：冷冻（-20℃）避光

宿 主：Rabbit

反应种属：Human Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the intercellular adhesion molecule (ICAM) family. All ICAM proteins are type I transmembrane glycoproteins, contain 2-9 immunoglobulin-like C2-type domains, and bind to the leukocyte adhesion LFA-1 protein. This protein may play a role in lymphocyte circulation by blocking LFA-1-dependent cell adhesion.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ICAM2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human ICAM2.
<b>Full name:</b>	intercellular adhesion molecule 2
<b>Synonyms :</b>	CD102

<b>SwissProt:</b>	P13598
<b>WB Predicted band size:</b>	31 kDa
<b>WB Positive control:</b>	NIH/3T3 cells lysate
<b>WB Recommended dilution:</b>	500-3000

