

FOXRED2 抗原（重组蛋白）

中文名称：FOXRED2 抗原（重组蛋白）

英文名称：FOXRED2 Antigen (Recombinant Protein)

别名：FAD dependent oxidoreductase domain containing 2; ERFAD

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 485-684 amino acids of human FOXRED2

技术规格

Full name:	FAD dependent oxidoreductase domain containing 2
Synonyms:	ERFAD
Swissprot:	Q8IWF2
Gene Accession:	BC040351
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Probable flavoprotein which may function in endoplasmic reticulum associated degradation (ERAD). May bind non-native proteins in the endoplasmic reticulum and target them to the ubiquitination machinery for subsequent degradation.