

## FTSJ3 抗原（重组蛋白）

中文名称：FTSJ3 抗原（重组蛋白）

英文名称：FTSJ3 Antigen (Recombinant Protein)

别名：FtsJ homolog 3; SPB1; EPCS3

储存：冷冻（-20℃）

相关类别：抗原

### 概述

Fusion protein corresponding to a region derived from 648-847 amino acids of human FTSJ3

### 技术规格

<b>Full name:</b>	FtsJ homolog 3
<b>Synonyms:</b>	SPB1; EPCS3
<b>Swissprot:</b>	Q8IY81
<b>Gene Accession:</b>	BC000131
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system</b> :	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli FtsJ protein and to a mouse protein possibly involved in embryogenesis.