

## FXR2 抗原（重组蛋白）

中文名称：FXR2 抗原（重组蛋白）

英文名称：FXR2 Antigen (Recombinant Protein)

别名：FMR1 autosomal homolog 2; FXR2P; FMR1L2

储存：冷冻（-20℃）

相关类别：抗原

### 概述

Fusion protein corresponding to a region derived from 474-673 amino acids of human FXR2

### 技术规格

<b>Full name:</b>	FMR1 autosomal homolog 2
<b>Synonyms:</b>	FXR2P; FMR1L2
<b>Swissprot:</b>	P51116
<b>Gene Accession:</b>	BC020090
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a RNA binding protein containing two KH domains and one RCG box, which is similar to FMRP and FXR1. It associates with polyribosomes, predominantly with 60S large ribosomal subunits. This encoded protein may self-associate or interact with FMRP and FXR1. It may have a role in the development of fragile X mental retardation syndrome.