

## 兔抗 ICOS 多克隆抗体

- 中文名称： 兔抗 ICOS 多克隆抗体
- 英文名称： Anti-ICOS rabbit polyclonal antibody
- 别名： inducible T cell costimulator; AILIM; CD278; CVID1
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 抗原： ICOS
- 反应种属： Human
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the CD28 and CTLA-4 cell-surface receptor family. It forms homodimers and plays an important role in cell-cell signaling, immune responses, and regulation of cell proliferation. [provided by RefSeq, Jul 2008]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ICOS
<b>Immunogen:</b>	Fusion protein of human ICOS
<b>Full name:</b>	inducible T cell costimulator
<b>Synonyms:</b>	AILIM; CD278; CVID1
<b>SwissProt:</b>	Q9Y6W8
<b>ELISA Recommended dil</b>	5000-10000

<b>ution:</b>	
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	100-200

