

GATD3B 抗原（重组蛋白）

中文名称： GATD3B 抗原（重组蛋白）

英文名称： GATD3B Antigen (Recombinant Protein)

别名： glutamine amidotransferase like class 1 domain containing 3B; KNP-I

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 42-268 amino acids of human GATD3B

技术规格

Full name:	glutamine amidotransferase like class 1 domain containing 3B
Synonyms:	KNP-I
Swissprot:	A0A0B4J2D5
Gene Accession:	BC003587
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	May play a role in promoting mitochondrial enlargement in cone photoreceptor cells.

mlbio 梅联生物
Good elisakit producers