

兔抗 ID3 多克隆抗体

中文名称： 兔抗 ID3 多克隆抗体

英文名称： Anti-ID3 rabbit polyclonal antibody

别名： HEIR-1; bHLHb25

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ID3

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a helix-loop-helix (HLH) protein that can form heterodimers with other HLH proteins. However, the encoded protein lacks a basic DNA-binding domain and therefore inhibits the DNA binding of any HLH protein with which it interacts.
Applications:	ELISA, IHC
Name of antibody:	ID3
Immunogen:	Synthetic peptide of human ID3
Full name:	inhibitor of DNA binding 3, dominant negative helix-loop-helix protein
Synonyms :	HEIR-1; bHLHb25

SwissProt:	Q02535
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and Human brain
IHC Recommend dilution:	50-200

